

## Diverse Variety of Human Blood Products & Biological Materials

- Large & Small Collections
- Affordable Prices
- Short Lead Times
- Custom Specifications
- High-Quality & Consistency

Contact us:

(215) 795-2838

orders@lampire.com



**LAMPiRE**  
BIOLOGICAL SUPPLIES

Donor: Whole Blood  
Human K2 EDTA  
Catalog #: 7203707  
Collection Date: 04MAY2018  
Volume: 1 unit  
Lot Number: LS 24 01116  
Store: 2-8 C  
[www.lampire.com](http://www.lampire.com)



**LAMPiRE** BIOLOGICAL LABORATORIES  
P.O. BOX 270 • PIPERSVILLE, PA 18947 • TEL: 215-795-2838 • FAX: 215-795-0237

**CAUTION: FOR RESEARCH USE ONLY**

Human Whole Blood in K2 EDTA

GENDER: Not Specified	BLOOD TYPE: Not Specified
CATALOG#: 7203707	STORAGE: 2-8°C
UNIT #: LS 24 01116	VOLUME: 1 Unit
COLLECTION DATE: 04MAY2018	EXP. DATE: 04JUN2018

TEST RESULTS: The donor unit used for the preparation of this product was collected at FDA licensed centers and has been tested using FDA licensed tests and found Negative for HBsAg, HBV NAT, HIV 1/2 Antibody, HIV NAT, HCV Antibody, HCV NAT, RPR, WNV NAT, HTLV I/II (Leukocyte products only), and Zika virus RNA (using an investigational NAT). Donor has tested negative for antibodies to T. cruzi either on the current donation or at least one previous donation.

Not for Use in Products Subject to License under Section 351 of the Public Health Service Act. Not for Transfusion.  
Registration # 2550014

## Human Blood Products & Materials

LAMPiRE offers human material and blood for research and diagnostic purposes. Human blood and biological products are collected from healthy donors. In addition to standard quality control tests, all donors of blood and blood derivatives are tested for HBsAg, HIV, HCV, and syphilis.

- Whole blood is manually drawn from the donor
- Normal plasma is collected by automated techniques
- Defibrinated plasma is prepared by adding calcium chloride to thawed citrated plasma, resulting in a clot. This is followed by separation of the serum as needed
- Sera from anticoagulant-free whole blood is prepared by collecting whole blood without anticoagulant, allowing clot formation, then separating the resulting serum

### Options to Meet Your Specifications:

- Gender
- Age
- Drug Free
- Specific Blood Type
- Non-Smoker/Smoker
- Bleed and Ship or with QC Testing

NOTE: Lampire's Human Blood and Biological Products are not for use in products subject to license under Section 351 of the Public Health Service Act and are not for use in transfusion.

